



## MOUSE-SIZED STAINLESS STEEL RUNNING WHEEL

Our stainless steel 11.5 cm (4.5 inch) diameter running wheel enables you to monitor voluntary running wheel behavior 24-hours-a-day while obtaining information on circadian rhythms and activity. This running wheel is designed specifically to meet your needs to use standard cages and racks and collect wheel turn data 24/7. Wheel turns are readily monitored using the magnetic reed switch.

### SPECIAL FEATURES:

- Stainless steel prevents rusting
- Bronze bushing for durability and noise-free operation
- No tools required for assembly or disassembly
- Entire assembly is suitable for commercial cage washers
- Wheels are completely contained within the cage – no custom cages or racks required
- Is filter-top compatible

### DESIGNED TO:

- Take up less lab space
- Enable a low-cost wheel solution
- Provide enrichment
- Be compatible with a digital counter
- Be autoclavable
- Be fully compatible with our VitalView data acquisition system and ER4000 Transponders